



Data Sheet

GUINEA PIG IgG

ANTIBODY, POLYCLONAL

Catalog no.:	G 4017
Immunogen:	Purified guinea pig IgG
Host:	Sheep
Matrix:	Serum
Specificity:	The antibody reacts specifically with guinea pig IgG
Contents:	10 ml
Known applications:	RIA (second antibody, 1:30)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at:	- 20 °C
	Repeated thawing and freezing must be avoided
Last updated on:	4 May 2022

For research use only

Publishing research using G 4017? Please let us know so that we can cite your publication as a reference.

