

Data Sheet

HUMAN ATRIAL NATRIURETIC PEPTIDE [ANP]

(aa 1-28)

SYNTHETIC

Catalog No.:	AE1020AG.1 / AE1020AG.2
Name:	Human ANP (aa 1-28) (SLRRSSCFGGRMDRIGASGLGCNSFRY)
Synonyms:	Atrial natriuretic factor (ANF)
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA, ANP, PND GeneID: 4878
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20°C Repeated thawing and freezing must be avoided
References:	1. Meyer M, Armbruster FP, Maier I, Schwing J, Haller C, Missbichler A, Adermann K, Forssmann W-G, Woloszczuk W (1998). Urinary and plasma urodilatin measured by a direct RIA using a highly specific antiserum. <i>Clin Chem</i> 44 (12): 2524-2529.
Last updated on:	12 May 2022

For research use only

Publishing research-using AE1020AG? Please let us know so that we can cite your publication as a reference.