

Data Sheet

HUMAN pro-ATRIAL NATRIURETIC PEPTIDE [ANP] (aa 12-25) SYNTHETIC

Catalog no.:	AE6201AG.1 / AE6201AG.2
Name:	Human pro-ANP (aa 12-25) (LMDFK <u>N</u> LLDHLEEK <u>C</u>) <small>Underlined amino acids differ from the original sequence; they were added or substituted for technical reasons.</small>
Synonyms:	Cardiodilatin-related peptide
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA, ANP, PND GeneID: 4878
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	19 May 2022

For research use only

Publishing research using AE6201AG? Please let us know so that we can cite your publication as a reference.

