



Data Sheet

BOVINE GRANULOCYTE MACROPHAGE COLONY-STIMULATING FACTOR [GM-CSF]

ANTIBODY, POLYCLONAL

Catalog no.:	AS1021.1/AS1021.2
Immunogen:	Recombinant bovine Granulocyte Macrophage Colony-stimulating Factor
Synonyms:	GM-CSF, Colony-stimulating factor
Swiss-Prot No:	P11052
Gene Information:	Gene Name: CSF2 GeneID: 281095.
Host:	Rabbit
Matrix:	Serum
Specificity:	Bovine GM-CSF
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA, Immunohistochemistry This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	24 May 2022

For research use only

Publishing research using AS1021? Please let us know so that we can cite your publication as a reference.

