



Data Sheet

BOVINE INTERFERON γ (IFN γ)

ANTIBODY, POLYCLONAL

Catalog no.:	AS1019.1/AS1019.2
Immunogen:	Recombinant bovine Interferon γ
Synonyms:	BoIFNG, IFN-gamma
Swiss-Prot No:	P07353
Gene Information:	Gene Name: IFNG GenelD: 281237
Host:	Rabbit
Matrix:	Serum
Specificity:	Bovine Interferon γ
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.

Known applications: Western Blot¹

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

References: 1. Kliem H, Berisha B, Meyer HH, Schams D (2009). Regulatory changes of apoptotic factors in the bovine corpus luteum after induced luteolysis. *Mol Reprod Dev* **76**(3): 220-230.

Last updated on: 28 April 2022

For research use only

Publishing research using AS1019? Please let us know so that we can cite your publication as a reference.

