



Data Sheet

CHICKEN EGG LYSOZYME NATIVE

Catalog no.:	A 1221AG.1 / A 1221AG.2
Name:	Lysozyme
Synonyms:	Mucopetide-N-acetylmuramoyl-hydrolase
Source:	Purified from chicken egg whites
Appearance:	White powder
Contents:	10 µg / 100 µg (lyophilized) Resuspend in aqua bidest. for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	04 May 2022

For research use only

Publishing research using A 1221AG? Please let us know so that we can cite your publication as a reference.

