

Data Sheet

BOVINE THYMOSIN β 9 (aa 33-41)

SYNTHETIC

Catalog no.:	A 9583AG.1 / A 9583AG.2
Name:	Bovine Thymosin β 9 (aa 33-41)
Sequence:	TIEQEKQAK
Synonyms:	Thymosin β -10
Swiss-Prot No:	P21752
Gene Information:	Gene Name: TMSB10 Gene ID: 282385
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g/1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9583AG? Please let us know so that we can cite your publication as a reference.

