



## Data Sheet

# BOVINE TNF $\alpha$

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AS1025.1 / AS1025.2
<b>Immunogen:</b>	Recombinant bovine Tumor Necrosis Factor alpha (TNF $\alpha$ )
<b>Synonyms:</b>	Tumor necrosis factor ligand superfamily member 2, Cachectin
<b>Swiss-Prot No:</b>	Q06599
<b>Gene Information:</b>	Gene Name: TNF, TNFA, TNFSF2 Gene ID: 280943
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Bovine TNF $\alpha$ No cross reaction with other growth factors
<b>Contents:</b>	20 $\mu$ l / 100 $\mu$ l (lyophilized) Resuspend in 20 $\mu$ l / 100 $\mu$ l aqua bidest.
<b>Known applications:</b>	RIA, Western Blot, immunohistochemistry  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 $^{\circ}$ C (lyophilized); - 20 $^{\circ}$ C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	28 April 2022

### For research use only

**Publishing research using AS1025? Please let us know so that we can cite your publication as a reference.**

