



Data Sheet

ANTI HUMAN TUMOR NECROSIS FACTOR- α [TNF α]

ANTIBODY, POLYCLONAL

Catalog no.:	AS1017.1 / AS1017.2
Immunogen:	Recombinant human tumor necrosis factor α (TNF α)
Reference No:	LocusID 7124
Host:	Rabbit
Matrix:	Serum
Specificity:	Human tumor necrosis factor α
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	RIA, Western Blot, immunohistochemistry This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 $^{\circ}$ C (lyophilized); - 20 $^{\circ}$ C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	28 April 2022

For research use only

Publishing research using AS1017? Please let us know so that we can cite your publication as a reference.

