


Data Sheet

1,25-DIHYDROXYVITAMIN D₃ LC-BIOTIN

SYNTHETIC

Catalog no.: A 1024AG.5 / A 1024AG.6
Name: Synthetic 1,25-dihydroxyvitamin D₃ LC-Biotin

Structure:



Mol. formula: C₄₆H₇₆N₄O₆S
Appearance: colourless, amorphous solid
Purity: >99 %
Solubility: Soluble in EtOH/MeOH and DMF, well soluble in halogenated hydrocarbons, not soluble in water
Contents: 100 ng / 1 µg (lyophilized)
Resuspend in 99% EtOH p.a. for required concentration
Store at: - 20 °C, protect from heat and light
References: Patent DE 198 40435 A1 registered on 25.06.1998
Last updated on: 04 May 2022

For research use only

Publishing research using A 1024AG? Please let us know so that we can cite your publication as a reference.

