



Data Sheet

25-HYDROXYVITAMIN D₃ LC-BIOTIN

SYNTHETIC

Catalog no.:	A 1025AG.5 / A 1025AG.6
Name:	Synthetic 25-hydroxyvitamin D ₃
Mol. formula:	C ₂₇ H ₄₄ O ₂
Synonyms:	Calcidiol; 25-hydroxycholecalciferol
Appearance:	Almost white crystals
Purity:	>99 %
Solubility:	Soluble in EtOH/MeOH and DMF, well soluble in halogenated hydrocarbons, not soluble in water
Contents:	100 ng / 1 µg (lyophilized) Resuspend in 99% EtOH p.a. for required concentration
Store at:	- 20 °C, protect from heat and light
Last updated on:	04 May 2022

For research use only

Publishing research using A 1025AG? Please let us know so that we can cite your publication as a reference.

