



Data Sheet

1,25-DIHYDROXYVITAMIN D₃

SYNTHETIC

Catalog no.:	A 1024AG.1 / A 1024AG.2
Name:	Synthetic 1,25-dihydroxyvitamin D ₃
Mol. formula:	C ₂₇ H ₄₄ O ₃
Synonyms:	Calcitriol, 1,25-dihydroxycholecalciferol
Appearance:	Almost white crystals
Purity:	>99 %
Solubility:	Soluble in EtOH/MeOH and DMF, well soluble in halogenated hydrocarbons, not soluble in water
Contents:	1 µg / 10 µg (lyophilized) Resuspend in 99% EtOH p.a. for required concentration
Store at:	- 20 °C, protect from heat and light
Last updated on:	04 May 2022

For research use only

Publishing research using A 1024AG? Please let us know so that we can cite your publication as a reference.

