



Data Sheet

HUMAN AMACO (aa 496-732)

ANTIBODY, POLYCLONAL

Catalog no.:	AP1013.1 / AP1013.2
Immunogen:	Recombinant human AMACO (aa 496-732), C-terminal His-tag
Synonyms:	Von Willebrand factor A domain-containing protein 2, A domain-containing protein similar to matrilin and collagen (AMACO), Colon cancer secreted protein 2 (CCSP-2)
Swiss-Prot No:	Q5GFL6
Gene Information:	Gene Name: VWA2, AMACO GenelD: 340706
Host:	Rabbit
Matrix:	Serum
Specificity:	Recombinant human Amaco (aa 496-732)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA, Western Blot, immunohistochemistry (paraffin sections) <small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.</small>
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	27 April 2022

For research use only

Publishing research using AP1013? Please let us know so that we can cite your publication as a reference.

