



## Data Sheet

# HUMAN IgG

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	G 4001
<b>Immunogen:</b>	Purified human IgG
<b>Host:</b>	Donkey
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human IgG
<b>Contents:</b>	10 ml
<b>Known applications:</b>	RIA (second antibody)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

<b>Store at:</b>	- 20 °C
	Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	4 May 2022

### For research use only

**Publishing research using G 4001? Please let us know so that we can cite your publication as a reference.**

