



## Data Sheet

# GOAT IgG

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	G 4002
<b>Immunogen:</b>	Purified goat IgG
<b>Host:</b>	Donkey
<b>Matrix:</b>	Serum
<b>Specificity:</b>	The antibody reacts specifically with goat IgG
<b>Contents:</b>	10 ml
<b>Known applications:</b>	RIA (second antibody)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

<b>Store at:</b>	- 20 °C
	Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	4 May 2022

### For research use only

**Publishing research using G 4002? Please let us know so that we can cite your publication as a reference.**

