



Data Sheet

HUMAN IgG

ANTIBODY, POLYCLONAL

Catalog no.:	G 4001
Immunogen:	Purified human IgG
Host:	Donkey
Matrix:	Serum, 137 mM D-(-)-mannitol
Specificity:	Human IgG
Contents:	10 ml (lyophilized) Resuspend in aqua bidest for required concentration.

Known applications: RIA (second antibody)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

Last updated on: 27 March 2025

For research use only

Publishing research using G 4001? Please let us know so that we can cite your publication as a reference.

