



Data Sheet

RABBIT IgG

ANTIBODY, POLYCLONAL

Catalog no.:	G 4004
Immunogen:	Purified rabbit IgG
Host:	Donkey
Matrix:	Serum, 137 mM D-(-)-mannitol
Specificity:	The antibody reacts specifically with rabbit IgG
Contents:	10 ml (lyophilized) Resuspend in aqua bidest for required concentration.
Known applications:	RIA (second antibody, 1:8) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Richter R, Bistran R, Escher S, Forssmann W-G, Vakili J, Henschler R, Spodsberg N, Frimpong-Boateng A, Forssmann U (2005). Quantum Proteolytic Activation of Chemokine CCL15 by Neutrophil Granulocytes Modulates Mononuclear Cell Adhesiveness. <i>The Journal of Immunology</i> 175 (3): 1599-1608.
Last updated on:	27 March 2025

For research use only

Publishing research using G 4004? Please let us know so that we can cite your publication as a reference.

