



Data Sheet

CHICKEN IgY

ANTIBODY, POLYCLONAL

Catalog no.:	G 4014
Immunogen:	Purified IgY
Host:	Donkey
Matrix:	Serum, 137 mM D-(-)-mannitol
Specificity:	The antibody reacts specifically with IgY
Contents:	10 ml (lyophilized) Resuspend in aqua bidest for required concentration.

Known applications: RIA (second antibody, 1:32)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

References: 1. Karmatschek M, Maier I, Seibel MJ, Woitge HW, Ziegler R, Armbruster FP (1997). Improved purification of human bone sialoprotein and development of a homologous radioimmunoassay. *Clin Chem* **43**(11): 2076-2082.

Last updated on: 27 March 2025

For research use only

Publishing research using G 4014? Please let us know so that we can cite your publication as a reference.

