



Data Sheet

APHRODISIN PEPTIDE (aa 107-126) MESOCRICETUS AURATUS (GOLDEN HAMSTER)

ANTIBODY, POLYCLONAL

Catalog no.:	AE1062.1 / AE1062.2
Immunogen:	Synthetic Aphrodisin Peptide (aa 107-126), poly-Lysine conjugated (KIFFTNKNMDRAGQETNMIV)
Swiss-Prot No:	Q9Z1I7
Host:	Rabbit
Matrix:	Serum
Specificity:	Aphrodisin Peptide (aa 107-126)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	Western Blot (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE1062? Please let us know so that we can cite your publication as a reference.

