



Data Sheet

HUMAN PRO-ATRIAL NATRIURETIC PEPTIDE [PRO-ANP] (aa 12-25)

ANTIBODY, MONOCLONAL

Catalog no.:	AE6204.1 / AE6204.2
Immunogen:	Synthetic human pro-ANP (aa 12-25) poly-Lysine conjugated (LMDFKNLLDHLEEK(C))
Synonyms:	Cardiodilatin-related peptide
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA GenelD: 4878
Host:	Mouse Balb/c
Clone no.:	A136-4-L1-23-G7
Isotype:	IgG ₁
Matrix:	Cell culture supernatant, Prot. G purified; PBS pH 7,4
Specificity:	Synthetic human pro-ANP (aa 12-25)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE6204? Please let us know so that we can cite your publication as a reference.

