

## Data Sheet

# HUMAN $\beta$ -DEFENSIN 3 (aa 6-22)

## SYNTHETIC

<b>Catalog no.:</b>	AE1050AG.1 / AE1050AG.2
<b>Name:</b>	Human $\beta$ -Defensin 3 (aa 6-22)
<b>Sequence:</b>	LQKYYCRVRGGRCVLS
<b>Synonyms:</b>	Beta-defensin 103, DEFB-3
<b>Swiss-Prot No:</b>	P81534
<b>Gene Information:</b>	Gene Name: DEFB103A and DEFB103B GeneID: 55894 and 4143255
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 $\mu$ g/ 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	18 May 2022

**For research use only**

**Publishing research using AE1050AG? Please let us know so that we can cite your publication as a reference.**

