

Data Sheet

HUMAN β -DEFENSIN 1 (aa 1-5)

SYNTHETIC

Catalog no.:	AE1047AG.1 / AE1047AG.2
Name:	Human β -Defensin 1 (aa 1-5)
Sequence:	DHYNC(ACM)GG
Synonyms:	β -Defensin 1 (BD-1), Defensin β 1
Swiss-Prot No:	P60022
Gene Information:	Gene Name: DEFB1 GeneID: 1672
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	18 May 2022

For research use only

Publishing research using AE1047AG? Please let us know so that we can cite your publication as a reference.

