



## Data Sheet

# HUMAN BONE SIALOPROTEIN [BSP] (aa 201-211)

## ANTIBODY, POLYCLONAL

**Catalog no.:** A 4223.1 / A 4223.2  
**Immunogen:** Synthetic human Bone Sialoprotein (aa 201-211) ((Y)TGGQGKGTSKT)  
**Synonyms:** Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein

**Swiss-Prot No:** P21815

**Gene Information:** Gene Name: IBSP  
 Gene ID: 3381

**Host:** Rabbit

**Matrix:** Serum

**Specificity:** Human Bone Sialoprotein (aa 201-211)  
 Cross reactivity with human Bone Sialoprotein

**Contents:** 20 µl / 100 µl (lyophilized)  
 Resuspend in 20 µl / 100 µl aqua bidest.

**Known applications:** ELISA (1:1000),  
 Western Blot (1:7000),  
 Immunocytochemistry (1:100)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

**Store at:** 2-8 °C (lyophilized); -20 °C (dissolved)  
 Repeated thawing and freezing must be avoided

**Last updated on:** 06 April 2022

**For research use only**

**Publishing research using A 4223? Please let us know so that we can cite your publication as a reference.**

