



Data Sheet

HUMAN BONE SIALOPROTEIN (BSP)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | A 1110.1 / A 1110.2 |
| Immunogen: | Recombinant human Bone Sialoprotein |
| Synonyms: | Bone sialoprotein 2, BSP2, Cell-binding sialoprotein, Integrin-binding sialoprotein |
| Swiss-Prot No: | P21815 |
| Gene Information: | Gene name: IBSP GeneID: 3381 |
| Host: | Goat |
| Matrix: | Serum |
| Specificity: | Human Bone Sialoprotein |
| Contents: | 20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest. |
| Known applications: | ELISA, RIA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 16 March 2022 |

For research use only

Publishing research using A 1110? Please let us know so that we can cite your publication as a reference.