

## Data Sheet

# HUMAN BRAIN NATRIURETIC PEPTIDE [BNP]

(aa 1-32)

SYNTHETIC

<b>Catalog no.:</b>	A 1015AG.1 / A 1015AG.2
<b>Name:</b>	Human BNP (aa 1-32), cyclic
<b>Synonyms:</b>	Natriuretic peptides B, Gamma-brain natriuretic peptide
<b>Sequence:</b>	SPKMVQSGSGCFGRKMDRISSSSSGLGCKVLRRH
<b>MW:</b>	3.5 kDa
<b>Swiss-Prot No:</b>	P16860
<b>Gene Information:</b>	Gene Name: NPPB GeneID: 4879
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg /1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	04 May 2022

**For research use only**

**Publishing research using A 1015AG? Please let us know so that we can cite your publication as a reference.**

