



Data Sheet

HUMAN C-PEPTIDE

ANTIBODY, POLYCLONAL, AFFINITY PURIFIED

Catalog no.:	A 5512.1 / A 5512.2
Immunogen:	C-Peptide conjugate EAEDLQVGQVKKKC-KLH
Synonyms:	Insulin
Swiss-Prot No:	P01308
Gene Information:	Gene Name: INS GenelD: 3630
Host:	Sheep
Preparation:	Affinity purified from sheep serum against target peptide and cross-adsorbed on an affinity column of immobilised mouse immunoglobulin
Matrix:	PBS pH7.4, 0.1% NaN ₃
Specificity	Human C-Peptide
Contents:	10µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	6 April 2022

For research use only

Publishing research using A 5512? Please let us know so that we can cite your publication as a reference.

