

Data Sheet

HUMAN PRO-C-TYPE NATRIURETIC PEPTIDE [proCNP] (aa 30-50) SYNTHETIC

Catalog No.:	AE1001AG.1 / AE1001AG.2
Name:	Human proCNP (aa 30-50)
Sequence:	GDKAPGGGGANLKGDRSRLLR
Swiss-Prot No:	P23582
Gene Information:	Gene Name: NPPC, CNP2 GeneID: 4880
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20°C Repeated thawing and freezing must be avoided
Last updated on:	12 May 2022

For research use only

Publishing research using AE1001AG? Please let us know so that we can cite your publication as a reference.

