



Data Sheet

HUMAN CALCITONIN GENE-RELATED PEPTIDE 2 SYNTHETIC

Catalog No.:	A 4612AG.1 / A 4612AG.2
Name:	Synthetic human Calcitonin gene-related peptide 2
Synonyms:	Beta-type CGRP, Beta-CGRP, Calcitonin gene-related peptide II, CGRP-II
Sequence:	ACNTATCVTHRLAGLLSRSGGMVKSNFVPTNVGSKAF
MW:	3.8 kDa
Swiss-Prot No:	P10092
Gene Information:	Gene Name: CALCB, CALC2 Gene ID: 797
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest. for required concentration.
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	05 May 2022

For research use only

Publishing research using A 4612AG? Please let us know so that we can cite your publication as a reference.

