



Data Sheet

HUMAN MRP8/14 (CALPROTECTIN, L1 PROTEIN)

ANTIBODY, POLYCLONAL

Catalog no.:	A 1075.1 / A 1075.2
Immunogen:	Recombinant human Calprotectin, tetrameric heterocomplex of MRP8 and MRP14
Synonyms:	Calprotectin, S100A8/A9, Calgranulin A/B heterocomplex, L1 Protein
Swiss-Prot No:	P05109 and P06702
Gene Information:	Gene Name: S100A8, MRP8 and S100A9, MRP14 GenelID: 6279 and 6280
Host:	Goat
Matrix:	Serum
Specificity:	Human calprotectin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	16 March 2022

For research use only

Publishing research using A 1075? Please let us know so that we can cite your publication as a reference.

