



Data Sheet

HORSE CHORIONIC GONADOTROPIN (CG)

ANTIBODY, POLYCLONAL

Catalog no.:	AS1004.1 / AS1004.2
Immunogen:	Purified horse Chorionic Gonadotropin
Synonyms:	Pregnant Mare Serum Gonadotropin (PMSG)
Host:	Rabbit
Matrix:	Serum
Specificity:	Horse Chorionic Gonadotropin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:500.000), Western Blot (1:500-10.000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions/concentrations in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	28 April 2022

For research use only

Publishing research using AS1004? Please let us know so that we can cite your publication as a reference.

