



## Data Sheet

# CREATINE KINASE (*SUBERITES DOMUNCULA*)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AC1009.1 / AC1009.2
<b>Immunogen:</b>	Purified Creatine kinase from <i>Suberites domuncula</i>
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	<i>Suberites domuncula</i> and <i>Geodia cydonium</i> Creatine kinase
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:1000), Western Blot (1:200), immunohistochemistry (1:20) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	14 April 2022

**For research use only**

**Publishing research using AC1009? Please let us know so that we can cite your publication as a reference.**

