



Data Sheet

CYCLIC GMP [cGMP]

ANTIBODY, POLYCLONAL

Catalog no.:	AE1132.1 / AE1132.2
Immunogen:	Synthetic cGMP, poly Lysine conjugated
Synonyms:	Cyclisches Guanosinmonophosphat
Host:	Rabbit
Matrix:	Serum
Specificity:	cGMP, no cross-reactivity with cAMP
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:800) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	22 April 2022

For research use only

Publishing research using AE1132? Please let us know so that we can cite your publication as a reference.

