



## Data Sheet

# HUMAN DERMCIDIN PEPTIDE (aa 96-110)

## ANTIBODY, POLYCLONAL

|                            |   |
|----------------------------|---|
| <b>Catalog no.:</b>        | AE1061.1 / AE1061.2   |
| <b>Immunogen:</b>          | Synthetic human Dermcidin Peptide (aa 96-110), polyLysine-conjugated (KGAVHDVKDVLDSVL)  |
| <b>Synonyms:</b>           | Preproteolysin  |
| <b>Swiss-Prot No:</b>      | P81605  |
| <b>Gene Information:</b>   | Gene Name: DCD<br>GeneID: 117159  |
| <b>Host:</b>               | Rabbit  |
| <b>Matrix:</b>             | Serum   |
| <b>Specificity:</b>        | Human Dermcidin Peptide (aa 96-110)   |
| <b>Contents:</b>           | 20 µl / 100 µl (lyophilized)<br>Resuspend in 20 µl / 100 µl aqua bidest.  |
| <b>Known applications:</b> | ELISA (1: 60 000), Western Blot (1:10 000)<br><small>This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.</small> |
| <b>Store at:</b>           | 2-8 °C (lyophilized); - 20 °C (dissolved)<br>Repeated thawing and freezing must be avoided  |
| <b>Last updated on:</b>    | 22 April 2022   |

### For research use only

**Publishing research using AE1061? Please let us know so that we can cite your publication as a reference.**

