



## Data Sheet

# PORCINE DIAMINE OXIDASE (DAO)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 1026.1 / A 1026.2
<b>Immunogen:</b>	Porcine Diamine oxidase
<b>Synonyms:</b>	Amiloride-sensitive amine oxidase, Amiloride-binding protein (ABP), Histaminase, PK-DAO
<b>Swiss-Prot No:</b>	Q9TRC7
<b>Gene Information:</b>	Gene Name: AOC1, ABP1 GeneID: 100517436
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Porcine DAO
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA, Western Blot (1:5000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	16 March 2022

### For research use only

**Publishing research using A 1026? Please let us know so that we can cite your publication as a reference.**

