



## Data Sheet

# HUMAN DIAMINE OXIDASE (DAO)

## ANTIBODY, MONOCLONAL

<b>Catalog no.:</b>	A 7004.1 / A 7004.2
<b>Immunogen:</b>	Recombinant human Diamine oxidase (DAO)
<b>Synonyms:</b>	Amiloride-sensitive amine oxidase (copper-containing), Amiloride-binding protein (ABP), Histaminase Kidney amine oxidase (KAO), EC 1.4.3.22
<b>Swiss-Prot No:</b>	P19801
<b>Gene Information:</b>	Gene Name: ABP1, AOC1, DAO1 GeneID: 26
<b>Host:</b>	Mouse Balb/c
<b>Clone no.:</b>	FP100
<b>Isotype:</b>	IgG <sub>1</sub>
<b>Matrix:</b>	Cell culture supernatant, Protein G purified, PBS pH 7.4
<b>Specificity:</b>	Diamine oxidase
<b>Contents:</b>	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved)  Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	7 April 2022

### For research use only

**Publishing research using A 7004? Please let us know so that we can cite your publication as a reference.**

