



Data Sheet

BOVINE ENDOTHELIN-1

ANTIBODY, POLYCLONAL

Catalog no.: AS1026.1 / AS1026.2

Immunogen: Bovine Endothelin-1

Synonyms: ET-1

Swiss-Prot No: P17322

Gene Information: Gene name: EDN1
GeneID: 281137

Host: Rabbit

Matrix: Serum

Specificity: Bovine Endothelin-1, human Endothelin-1

Contents: 20µl / 100 µl (lyophilized)
Resuspend in aqua bidest. for required concentration

Known applications: RIA, ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)
Repeated thawing and freezing must be avoided

Last updated on: 15 November 2024

For research use only

Publishing research using AS1026? Please let us know so that we can cite your publication as a reference.

