

Data Sheet

HUMAN EPIDERMAL GROWTH FACTOR RECEPTOR (EGF-R)

NATIVE

Catalog no.:	A 1016AG.1 / A 1016AG.2
Name:	Human EGF-Receptor
Synonyms:	Receptor tyrosine-protein kinase erbB-1, Proto-oncogene c-ErbB-1
Swiss-Prot No:	P00533
Gene Information:	Gene Name: EGFR, ERBB, ERBB1, HER1 Gene ID: 1956
Source:	Human epidermal cell line
Purification:	Affinity chromatography
Storage buffer:	PBS, 50 µg BSA / 1 µg cytokin
Purity:	> 95% (SDS-PAGE)
Contents:	10 µg / 100 µg (lyophilized) Resuspend in PBS / 0.1% BSA for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	4 May 2022

For research use only

Publishing research using A 1016AG? Please let us know so that we can cite your publication as a reference.