

## Data Sheet

# HUMAN GASTRIN

## SYNTHETIC

<b>Catalog No.:</b>	A 5511AG.1 / A 5511AG.2
<b>Name:</b>	Human Gastrin
<b>Synonyms:</b>	Gastrin component III, Gastrin-17
<b>Sequence:</b>	Pyr-GPWLEEEEEAYGWMDF ( <i>Pyr= Pyroglutamic Acid</i> )
<b>MW:</b>	2.12 kDa
<b>Swiss-Prot No:</b>	P01350
<b>Gene Information:</b>	Gene Name: GAST, GAS Gene ID: 2520
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 60 % acetonitrile in 0.01 M ammonium bicarbonate buffer for required concentration.
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	6 May 2022

### For research use only

**Publishing research using A 5511AG? Please let us know so that we can cite your publication as a reference.**

