



Data Sheet

HUMAN GLUCAGON-LIKE PEPTIDE-2 [GLP-2]

ANTIBODY, POLYCLONAL

Catalog no.:	A 6102.1 / A 6102.2
Immunogen:	Synthetic human GLP-2, KLH-conjugated (HADGSFSDEMNTILDNLAARDFINWLIQTKITD)
Synonyms:	GLP-2, Glucagon (aa 126–158)
Swiss-Prot No:	P01275
Gene Information:	Gene Name: GCG GenelD: 2641
Host:	Rabbit
Matrix:	Serum
Specificity:	Human GLP-2
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:3000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	6 April 2022

For research use only

Publishing research using A 6102? Please let us know so that we can cite your publication as a reference.

