



Data Sheet

HUMAN GLUCAGON-LIKE PEPTIDE-1 (GLP-1) SYNTHETIC

Catalog No.:	A 6104AG.1 / A 6104AG.2
Name:	Human GLP-1 Peptide (aa 7-37)
Sequence:	HAEGTFTSDVSSYLEGQAAKEFIAWLVKGRG
M_r:	3355.71
Swiss-Prot No:	P01275
Gene Information:	Gene Name: GCG GeneID: 2641
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	06 May 2022

For research use only

Publishing research using A 6104AG? Please let us know so that we can cite your publication as a reference.

