



Data Sheet

GLYPICAN-4 (aa 367-382)

ANTIBODY, POLYCLONAL

Catalog no.:	AH1009.1 / AH1009.2
Immunogen:	Synthetic glypican-4 (aa 367-382) coupled to KLH (RPHHPEERPTTAAGTSL)
Synonyms:	K-glypican
Swiss-Prot No:	O75487
Gene Information:	Gene Name: GPC4 GenelD: 2239
Host:	Rabbit
Matrix:	Serum
Specificity:	Specific for a core protein epitope (aa 367-382) of human glypican-4. Due to the highly conserved structure and sequence of glypican-4, it might react also with glypican-4 from other species, like mouse.
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	Western Blot (1:1000), immunohistochemistry (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Siebertz B, Drzeniek Z, Stöcker G, Just U, Handt S, Haubeck HD (1999): Expression of glypican-4 in haematopoietic progenitor and bone marrow stromal cells. <i>Biochem J</i> 344 : 937-4.
Last updated on:	27 April 2022

For research use only

Publishing research using AH1009? Please let us know so that we can cite your publication as a reference.

