



Data Sheet

GLYPICAN-4 (aa 509-520)

ANTIBODY, POLYCLONAL

Catalog no.:	AH1010.1 / AH1010.2
Immunogen:	Synthetic glypican-4 (aa 509-520) coupled to KLH (FEFDYNATDHAG)
Synonyms:	K-glypican
Swiss-Prot No:	O75487
Gene Information:	Gene Name: GPC4 GeneID: 2239
Host:	Rabbit
Matrix:	Serum
Specificity:	Human WB, IHC-p, mouse glypican-4
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	Western Blot (1:1000) ¹ , flowcytometry ¹ , immunohistochemistry (paraffin sections, 1:1000) ¹ This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Siebertz B, Drzeniek Z, Stöcker G, Just U, Handt S, Haubeck HD (1999). Expression of glypican-4 in haematopoietic progenitor and bone marrow stromal cells. <i>Biochem J</i> 344 : 937-43.
Last updated on:	27 April 2022



For research use only

Publishing research using AH1010? Please let us know so that we can cite your publication as a reference.

