

## Data Sheet

# HUMAN GUANYLIN (aa 80-94)

## SYNTHETIC

<b>Catalog No.:</b>	AE1017AG.1 / AE1017AG.2
<b>Name:</b>	Human Guanylin (aa 80-94) active peptide (PGTCEICAYAACTGC)
<b>Synonyms:</b>	Guanylate cyclase activator 2A, Guanylate cyclase-activating protein 1, Guanylate cyclase-activating protein I (GCAP-I)
<b>Swiss-Prot No:</b>	Q02747
<b>Gene Information:</b>	Gene Name: GUCA2A, GUCA2 GenelD: 2980
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	12 May 2022

**For research use only**

**Publishing research using AE1017AG? Please let us know so that we can cite your publication as a reference.**