


Data Sheet

INDOLEAMINE 2,3-DIOXYGENASE 1 (IDO-1)

ANTIBODY, POLYCLONAL

Catalog no.:	A 2106.1 / A 2106.2
Immunogen:	Recombinant human Indoleamine 2,3-Dioxygenase 1
Synonyms:	Indoleamine-pyrrole 2,3-dioxygenase, EC 1.13.11.52
Swiss-Prot No:	P14902
Gene Information:	Gene Name: IDO1, IDO, INDO GenelD: 3620
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Indoleamine 2,3-Dioxygenase 1
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:250) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	18 March 2022



For research use only

Publishing research using A 2106? Please let us know so that we can cite your publication as a reference.

