

## Data Sheet

# HUMAN KATACALCIN SYNTHETIC

<b>Catalog No.:</b>	A 1211AG.1 / A 1211AG.2
<b>Name:</b>	Human Katalcalcin (PDN-21)
<b>Synonyms:</b>	Calcitonin carboxyl-terminal peptide (CCP), PDN-21
<b>Sequence:</b>	DMSSDLERDHRPHVSMQPQAN
<b>Swiss-Prot No:</b>	P01258
<b>Gene Information:</b>	Gene Name: CALCA, CALC1 Gene ID: 796
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest for required concentration.
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	4 May 2022

**For research use only**

**Publishing research using A 1211AG? Please let us know so that we can cite your publication as a reference.**