



Data Sheet

HUMAN KATACALCIN

ANTIBODY, POLYCLONAL

Catalog no.:	A 1211.1 / A 1211.2
Immunogen:	Synthetic human Katalcalcin (PDN-21), Ovalbumin conjugated (DMSSDLERDHRPHVSMQAN)
Synonyms:	Calcitonin carboxyl-terminal peptide (CCP), PDN-21
Swiss-Prot No:	P01258
Gene Information:	Gene Name: CALCA, CALC1 GenelD: 796
Host:	Goat
Matrix:	Serum
Specificity:	Human Katalcalcin (PDN-21)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:500) ¹ This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Blind E, Raue F, Klaiber T, Zink A, Schroth J, Buhr H, Ziegler R (1992). Evaluation of sensitive PDN-21 (katalcalcin) determination as tumor marker in medullary thyroid carcinoma. <i>J Endocrinol Invest</i> 15 :93-8.
Last updated on:	16 March 2022

For research use only

Publishing research using A 1211? Please let us know so that we can cite your publication as a reference.

