

Data Sheet

HUMAN MATRIX Gla PROTEIN (MGP) SYNTHETIC

Catalog no.:	A 1009AG.1 / A 1009AG.2
Name:	human Matrix Gla Protein (aa 3-15)
Synonyms:	Cell growth-inhibiting gene 36 protein
Sequence:	SHESMESYELNPF
MW:	1569 Da
Swiss-Prot No:	P08493
Gene Information:	Gene Name: MGP, MGLAP Gene ID: 4256
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in aqua bidest for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	4 May 2022

For research use only

Publishing research using A 1009AG? Please let us know so that we can cite your publication as a reference.

