



## Data Sheet

# HUMAN MEPRIN A ALPHA (aa 649-654)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AE1136.1 / AE1136.2
<b>Immunogen:</b>	Synthetic human Meprin A alpha (aa 649-654) KLH conjugated (VHGSVL)
<b>Synonyms:</b>	Meprin A subunit alpha, Endopeptidase-2, N-benzoyl-L-tyrosyl-P-amino-benzoic acid hydrolase subunit alpha, PABA peptide hydrolase, PPH alpha
<b>Swiss-Prot No:</b>	Q16819
<b>Gene Information:</b>	Gene Name: MEP1A GenelD: 4224
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Synthetic human Meprin A alpha (aa 649-654)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	22 April 2022

### For research use only

**Publishing research using AE1136? Please let us know so that we can cite your publication as a reference.**

