



Data Sheet

HUMAN MYOSTATIN

ANTIBODY, POLYCLONAL

Catalog no.:	AK3009.1 / AK3009.2
Immunogen:	Recombinant human Myostatin
Synonyms:	Growth / differentiation factor 8 (GDF-8)
Swiss-Prot No:	O14793
Gene Information:	Gene Name: MSTN, GDF8 GeneID: 2660
Host:	Goat
Matrix:	Serum
Specificity:	Human Myostatin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	27 April 2022

For research use only

Publishing research using AK3009? Please let us know so that we can cite your publication as a reference.

