



Data Sheet

HUMAN PRO-MYOSTATIN (aa 79-92)

SYNTHETIC

Catalog no.:	AK 3004AG.1 / AK 3004AG.2
Name:	Human Pro-Myostatin (aa 79-92)
Sequence:	SDGSLEDDDYHATT
Synonyms:	Growth/differentiation factor 8 (GDF-8)
Swiss-Prot No:	O14793
Gene Information:	Gene name: MSTN, GDF8 GeneID: 2660
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	20 May 2022

For research use only

Publishing research using AK3004AG? Please let us know so that we can cite your publication as a reference.

